

# Disposable Specimen Collection Swab

## User Manual

Version: A/O

Edit Date: 2020-2-19

### 【Product Name】

Disposable Specimen Collection Swab

### 【Model Number】

DSC-KC、DSC-DKC、DSC-KE、DSC-DKE、TFS-TC、FS-H10、NFS-TB、NFS-TB1、TFS-TD、  
TFS-T、NFS-T、TFS-TE、TFS-TF、TFS-TG、UFS-TA、CFS-TH

### 【Description】

The swab is composed of the handle and the sampling head.

### 【Application】

Collect the microbiological samples from oral, nose, throat, cervix, urinary tract and other parts of human body.

### 【Contraindication】

Use with caution for those allergic to flocking.

### 【Product performance】

After irradiation sterilization, the product should be sterile.

### 【Instructions】

1. Open the package, take out the swab.
2. Sampling the specimen from human body.



**【Attentions】**

1. Make sure the good sealing of the package. If the package is damaged or has expired, it shall be prohibited to use.
2. Disposable use for one person.
3. Irradiation sterilized, valid for 2 years.

**【Storage】**

The store should be ventilated and at Room temperature.

No stress. No exposing to the sun or rain.

**【Manufacturer】** Shenzhen J.able Bio Co.,Ltd.

**【Address】** 4th floor, 13th Building, Lihao Industrial Park, Ainan Road #78, Longgang District, Shenzhen, Guangdong Province, China

**【Tel】** 0755-28591602

**【E-mail】** info@jablebio.com

**【Postal Code】** 518116